

Are you hearing a ringing or buzzing sound? Or maybe your hearing isn't what it used to be?



is now offering sound therapies developed by **AudioCardio**



AudioCardio™ Mobile App

AudioCardio™ is a mobile app that generates a personalized sound therapy (just inaudible) designed to protect and strengthen your hearing.

Threshold Sound Conditioning™

Threshold Sound Conditioning™ technology is backed by over 15 years of research and development by technologists, award winning music composers and scientists.

Backed By Data and Science

Threshold Sound Conditioning™ technology has been clinically tested and proven with Stanford University, Palo Alto Medical Foundation, Samsung Medical Center, Chung-Ang and Nanjing Universities.



Make it Part of Your Daily Routine

Use AudioCardio™ once a day for one hour while you enjoy other activities like walking, working out, listening to music or even while you work.

Real Results for Better Hearing Health

Some individuals reported noticeable results in as little as a week and over 70% heard significant changes after two weeks of consecutive use.

Subscribe and Save

AudioCardio™ is a subscription based service that can be purchased on a monthly, semi annual or annual subscription starting at \$8.33 USD per month.



Use promo code **EAR30**
at <https://app.audiocardio.com/signup/earinc>



You're not alone in handling this problem. Over 465 million people around the world experience some form of hearing loss or tinnitus. AudioCardio aims to provide an affordable and state of the art solution to managing this problem.

Best Practices

- 1 Connect your earbuds or headphones and go to a very quiet room.
***Avoid rooms with ambient noise (appliances, passing cars, etc.)**
- 2 Move the slider up until you can comfortably hear the tone being played and then move and set it at the point where you can barely hear it.
- 3 You will do this a total of ten times. First you will do this for the left ear, then the right.



- 4 After completing the assessment, select the "Start Therapy" button to start your personalized sound therapy. You may or may not hear faint frequencies coming from each side.
- 5 Enjoy your sound therapy for one hour per day while doing your favorite everyday tasks like cooking, watching TV, and exercising.

New product and service

by



and



Use promo code

EAR30

at <https://app.audiocardio.com/signup/earinc>

Learn more at www.audiocardio.com - Inquiries at info@audiocardio.com